

Silicone Tray Adhesive

Tray adhesives are recommended to assure reliable adhesion of impression material to the impression tray; and to avoid uncontrollable deformation during impression taking.



NOTE: The adhesive does not always secure sufficient adhesion to trays made of other plastic materials or individual trays made of base plates. A polyether adhesive must be used for polyether impression materials.

INSTRUCTIONS FOR USE

- 1. Apply a thin layer of Access® Tray Adhesive to the tray using the applicator brush.
- 2. Allow to dry for at least 3 minutes.
- 3. Load the tray with impression material and take the impression.

SAFETY DATA SHEET (MSDS)

In accordance with Hazard Communication Standard (HCS), 29CFR 1910.1200(g), 2012 revision

SECTION I - Product / Company Identification Product Trade Name: Access Tray Adhesive

Part (Item) Number: 360005

Description: Impression Tray Adhesive

Manufacturer: CENTRIX, INC. 770 River Road, Shelton, CT 06484 USA

Emergency Telephone Number: (203) 929 - 5582 Information Telephone Number: (800) 235 - 5862

SECTION II - Hazardous Ingredients/Identity Information

Designation of Risk: Xi, irritant Special Details of Risk (R-categories): 36-38

CLASSIFICATION: Flammable liquid. May cause dizziness, headache, nausea and mental confusion if inhaled. May cause irritation or dermatitison skin. May cause severe eye irritation. Harmful if ingested with both acute and chronicle affect.

SECTION III - Hazardous Ingredients

Adhesive for vinvl silicone

Hazardous Components: Alkyldisiloxane (107-45-0), Xylene (1130-20-7)

SECTION IV - First Aid Measures

Skin: Wipe off with towel. Wash skin with plenty of water and soap. If irritation develops, seek medical attention.

Ingestion: Seek medical attention immediately

Inhalation: Remove to fresh air. Seek medical attention if adverse health effects are experienced.

Eves: Flush immediately with plenty of water (for at least 15 minutes). If irritation persists, see eye doctor,

SECTION V - Fire-Fighting Measures

Suitable extinguishing media: Foam, powder, carbon dioxide

Forbidden extinguishing media: None

Special firefighting measures: None

Unusual fire and explosion hazards: Flammable material

Special protection equipment: Appropriate protective equipment must be worn in case of fire.

Solvent vapors may form explosive mixtures with air. Vapors may travel to source of ignition and flash back.

SECTION VI - Accidental Release Measures

Personal precautions: Eliminate all ignition sources. Contain and absorb spilled material using an appropriate solvent spill absorbent. Collect in a suitable container for disposal.

Environmental precautions: Residual material should evaporate from impervious surfaces. Do not flush in sewer system Reclaiming methods: Material can be burned in approved incinerator in accordance with federal, state and local regulations.

SECTION VII - Handling and Storage

Handling: Product should only be handled by dentists or upon their instructions. Wear appropriate protective equipment and clothing, Provide adequate ventilation, Prevent contact with skin and eyes, Avoid inhalation or ingestion,

Storage: Store in a cool safe area away from ignition sources. Keep container closed at all times. Keep away from heat, flame, or sparks.

Containers: Store in supplied packaging (keep container tightly closed).

Incompatibilities: None

Other precautions: None

SECTION VIII - Exposure Controls / Personal Protection

General measures: Ensure adequate ventilation when handling

Respiratory protection: Provide adequate ventilation.

Eve protection: Eve protection is recommended.

Skin protection: Wear appropriate protective equipment and clothing.

Hands protection: Impervious gloves are recommended.

Other equipment: None

SECTION IX - Physical and Chemical Properties

Appearance: Viscous fluid, red

Odor: Aromatic

Water solubility: Insoluble in water, miscible with most organic solvents

Boiling point: >80 C Flash point ca. -9 C

SECTION X - Stability and Reactivity

Stability: Stable

Conditions to avoid:

Materials to avoid: Avoid all potential ignition sources

Hazardous decomposition products: None to our knowledge Hazardous polymerization products: None to our knowledge

$C \in$

SECTION XI - Toxicological Information

Ingestion: Hazard Eve contact: Irritant Skin contact: Irritant Inhalation: Hazard

SECTION XII - Ecological Information

No data available

SECTION XIII - Disposal Considerations

Handle in accordance with applicable federal, state and local regulation. May be burned in approved incinerator.

SECTION XIV - Transport Information

This is a flammable material according to DOT regulations.

SECTION XV - Regulatory Information

 $Hazard\ Classification\ Symbol:\ This\ product\ contains\ the\ following\ toxic\ chemicals\ subject\ to\ reporting\ requirements:$ Xylene (CAS#1330-20-7)

Federal -TSCA: All components are either listed or otherwise comply with requirements.

International - DSL (Canada): All ingredients are listed on the Canadian Domestic

Substance List (DSL) International

WHMIS (Canada): Not a "controlled product" under the Canadian Workplace Hazardous Materials Information System- (WHMIS)

SECTION XVI - Other Information

Product handling exclusively by dental specialists per instructions for use.

CAUTION: PRODUCT FOR PROFESSIONAL USE ONLY.

The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.



Making Dentistry Easier.sm



770 River Road, Shelton, CT 06484 USA (203) 929-5582 (800) 235-5862 FAX (203) 929-6804 www.centrixdental.com



DENTAL-CONTACT VERTRIEBS OHG,

Im Spiet 101, DE 26506 Norden

CENTRIX INCORPORATED

EMERGO AUSTRALIA

Level 20, Tower II, Darling Park, 201 Sussex Street Sydney, NSW 2000, Australia

©2015 Centrix Inc. Made in U.S.A.

012022 REV D

Revised: 4/2015