Connexio[®]

Self-Etch Bond

English

INSTRUCTIONS FOR USE

For professional use only

Connexio Self-Etch Bond is a dual-curing, self-etching bonding agent reinforced with nano-particles for creating a durable bond between tooth substance and light-, self- and dual-curing restoratives without marginal leakage. This all-in-one adhesive consists of 2 components, which eliminates the need for a separate etching step, yet still achieves the same adhesive properties as total-etch bonds. Connexio Self-Etch Bond tolerates residual moisture and can be used with a variety of light, self and dual-curing composites (see contraindications).

INDICATIONS:

- Direct self- or dual-curing composite restorations and core buildups
- Direct light-curing composite-based restorations
- Fixing root pins/posts with dual- or self-curing composite cements
- Indirect restorations; with the application of dual-curing/self-curing composite cements for bonding inlays, onlays, crowns and bridges

CONTRAINDICATIONS:

- 1.** Connexio Self-Etch Bond contains organic acids, BIS-GMA, HEMA, TMPTMA, campherchinon, amines (DABE), BHT, catalysts, fluorides and ethanol. Do not use in patients with known allergies to these ingredients.
- 2.** Check instructions for use with any core buildup material for additional information.
- 3. Warnings and precautionary measures:**
 - Phenolic substances (e.g. temporary materials containing eugenol) disturb polymerization and therefore should not be used with Connexio Self-Etch Bond.
 - Avoid contact with the oral mucosa. Contact of the oral mucosa with Connexio may result in a whitish discoloration caused by the coagulation of proteins. This is a reversible occurrence that subsides within a few days. In case of contact with the eyes, rinse with plenty of clean water and consult a physician.
 - Application brushes provide precise and thin application of Connexio Self-Etch Bond.

INSTRUCTIONS FOR USE:

- 1. PREVENTING CONTAMINATION** – A dry working field is essential due to the subsequent application of the restorative. Avoid contamination with blood or saliva. Using a rubber dam is recommended.
- 2. PREPARATION** – Prepare the cavity according to the principles of the adhesive restoration technique. Rinse cavity thoroughly with water. Remove excess moisture with a faint air jet. Do not over dry dentin. The dentinal surface of the cavity should remain slightly moist. If the tooth has not been prepared (e.g. cervical defect) the surface must be thoroughly cleaned and enamel margin roughened with a finishing diamond. Areas in close proximity to the pulp should be protected with a suitable calcium hydroxide-based liner.
- 3. MIXING** – Connexio Self-Etch Bond is dual-curing and therefore must be applied immediately after mixing. Mix 1 drop of Part 1 and 1 drop of Part 2 on a mixing palette with an application brush for approximately 2 sec (produces a self-etching adhesive).
- 4. CONNEXIO SELF-ETCH BOND APPLICATION** – Apply the adhesive in a layer of medium thickness to the enamel/dentin and rub into the tooth surface for 20 sec.

CAUTION: The mixed bonding is light curing and thus an intensive exposure to ambient light should be avoided. The OP-Light during the application should be dimmed.

5. DRYING - Dry the adhesive layer for at least 5 sec with an air syringe.

6. POLYMERIZATION –

a. Direct Restoration

Light, self and dual-curing materials: The bonding layer must be polymerized with blue light (halogen or LED lights) for 10 sec when using light, self and dual-curing restoratives.

b. Indirect restoration

Light-curing luting composites: The bonding layer must be polymerized with blue light (halogen or LED lights) for 10 sec when using light-curing luting materials.

Self-curing or dual-curing luting materials: When using dual or self-curing luting materials, Connexio Self-Etch Bond also provides excellent adhesion without light-curing. Connexio Self-Etch Bond can also be light-cured for 10 sec with blue light (halogen or LED lights) as an optional step.

Luting root posts:

(a) Massage the mixed Connexio Self-Etch Bond (Part 1 and 2) into the occlusal surface of the core and root canal for 20 sec with an applicator. Disperse solvent with an oil-free air-jet for 5 sec. Remove excess material by dabbing it off with an absorbent strip of paper. Do not light-cure the bonding layer.

(b) Having prepared the root post in accordance with the manufacturer’s instructions, insert it into the root canal using a suitable luting material. The aim is to obtain a slight excess of cement.

(c) Light-cure for at least 40 sec to fix the root post and cure the Connexio Self-Etch Bond applied to the occlusal surface.

7. ADVICE ON LUTING MARYLAND BRIDGES – For luting Maryland bridges, the enamel must either be roughened or conditioned with a phosphoric acid etching gel before applying Connexio Self-Etch Bond.

Self-etching bonding materials are less effective on unprepared enamel. Applying bonding material to enamel that has not been ground or beveled can lead to margin discoloration. For optimum adhesion results with Connexio Self-Etch Bond, the enamel should be prepared in advance in class III, IV, V cavities and with diastema closures.

Advice for storage and usage:

Tightly close Connexio Self-Etch Bond bottles after use. Store in an upright position to guarantee flowback of the liquid. Store at temperatures between 39°F - 73°F (4°C - 23°C). Allow Connexio Self-Etch Bond to reach room temperature before use. Do not expose Connexio Self-Etch Bond to direct sunlight or the chair light.

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