# Champ®

## Pit & Fissure Sealant

Champ 16 x 0.25 mL Lollipop Pack, Champ 64 x 0.25 mL Lollipop Pack, and Champ 1 x 1.2 mL Syringe

## GB - ENGLISH

**INSTRUCTIONS FOR USE**

**CHAMP® PIT & FISSURE SEALANT**

Champ is intended to seal the pits and fissures in tooth enamel to aid in the prevention of caries.

**INSTRUCTIONS FOR USE:**

- Thoroughly clean enamel surfaces with oil-free cleaning paste. Remove all debris.
- Rinse well, removing all residual paste from pits and fissures in enamel.
- Isolate treatment area and dry lightly with oil-free compressed air.
- Apply phosphoric acid etching gel to the clean tooth surface for 15 seconds. Rinse well.
- Lightly dry and remove excess water with a cotton pellet or clean compressed air. Leave tooth surfaces slightly moist. Champ bonds to surfaces that are slightly moist; however, it is best to not let saliva contaminate the prep area.
- To open:
  - If using the single-dose Champ LollipopPack®, grip the foil cover and slowly peel back, exposing the brush applicator and sealant well.
  - If using the Champ syringe, remove cap, screw on applicator tip.
- Apply a thin film of Champ onto treatment area. Apply the material from cusp to cusp, but do not cover the marginal ridges.
- Light cure. Champ cures with all lights. Curing time for a halogen light with a minimum output intensity of 300 mW/cm2 is 20 seconds. More powerful lights will cure faster.
- Check and adjust occlusion if necessary.
- Throw away all packaging and single-dose LollipopPack once procedure complete.

**CONTRAINDICATIONS:**

Champ is contraindicated in patients who are allergic to any of its ingredients.

Do not use on dental caries.

Champ 16 x 0.25 mL Lollipop Pack, Champ 64 x 0.25 mL Lollipop Pack, and Champ 1 x 1.2 mL Syringe

Champ 16 x 0.25 mL Lollipop Pack, Champ 64 x 0.25 mL Lollipop Pack, and Champ 1 x 1.2 mL Syringe

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#### SAFETY DATA SHEET (SDS)

In accordance with Hazard Communication Standard (HCS), 29CFR 1910.1200(g), 2012 revision

#### SECTION I - PRODUCT / COMPANY IDENTIFICATION

REF - See Above

Company: Centrix Inc., 770 River Road, Shelton, CT 06484

Emergency / Information Telephone Numbers: 203-929-5582

Toll Free: 800-235-5862 (USA & Canada)

#### SECTION II – HAZARDOUS INGREDIENTS/IDENTITY INFORMATION

Classification: Irritant

**Adverse human health effects:** May be irritating to eyes, respiratory system and skin. May cause sensitization by skin contact.

Hazard Class	Hazard Category	Hazard Statement
Eye Irritation	2	H319
STOT SE	3	H335
Skin irritation	2	H315
Skin sensitization	1	H317

#### SECTION III - COMPOSITION

**Chemical characterization of the preparation:** Acrylate ester monomers in a glass-filled, light-cured paste. Hazardous ingredients

CAS Number	Name of the Ingredient	%	Classification per 67/548/EEC	Classification per Regulation (EC) No. 1278/2008 (CLP)
Proprietary R 3637//38 - 43	Uncured acrylate Eye Irritation, 2, H319	55-60%	Irritant; Xi	ester monomers <p>STOT SE 3, H335 Skin irritation, 2, H315</p>
112945-52-5	Silica, amorphous	5%	Irritant; Xi, R 36/37/38	Skin sensitization, 1, H317 <p>Eye irritation, 2, H319 STOT SE 3, H335 Skin irritation, 2, H315</p>

#### SECTION IV - FIRST AID MEASURES

**General Information:** May be irritating to eyes, respiratory system and skin. Do not inhale vapors. Avoid contact with eyes and skin. May cause sensitization by skin contact. Show this safety data sheet to medical personnel.

Get medical attention in case of uncertainty.

**Inhalation:** Move to fresh air. If necessary, administer oxygen and/or artificial respiration. Seek medical attention.
**Skin Contact:** Take off contaminated clothing. Wash skin thoroughly with soap and running water for 15 minutes.

**Eye Contact:** Keep eyelids apart and flush with running water for 15+ minutes. Get medical attention.

**Ingestion:** Rinse mouth with water. Do not induce vomiting. Get immediate medical attention. May be irritating to mucous membranes. Never give anything by mouth to an unconscious person.

**Precautions for first responders:** Ventilate the area. Wear eye and skin protection. Information for physicians
**Symptoms:** Irritation or redness in eyes, throat or on skin.

**Hazards:** May be irritating to eyes, respiratory system and skin. May cause sensitization by skin contact.

**Treatment:** As above under First Aid.

#### SECTION V - FIRE-FIGHTING MEASURES

**Suitable extinguishing media:** Carbon dioxide, dry chemical, alcohol foam, or water fog. Water spray may be used to keep fire exposed containers cool.

**Extinguishing media to avoid:** Do not use direct water stream

**Special exposure hazards in a fire:** Heat may cause polymerization with rapid release of energy.

**Special protective equipment for fire-fighters:** A self-contained breathing apparatus should be worn by fire fighting personnel. Cool exposed containers with water spray to prevent polymerization under fire conditions.

#### SECTION VI - ACCIDENTAL RELEASE MEASURES

**Personal precautions:** Wear safety glasses, gloves and lab coat.

**Environmental precautions:** Avoid releasing large amounts of uncured material into the environment. Cured / set-up material is, to our knowledge, inert.

**Method for clean up:** Contain spilled material. Absorb or wipe up spill with suitable material (paper towels or cloths). Collect for disposal in a covered container. Wash area of spill with alcohol or soap and water.

#### SECTION VII - HANDLING AND STORAGE

**Handling:** For use only by dental professionals. Follow good hygiene practices. Do not smoke, eat or drink while using. Avoid direct / strong light sources, extremes of temperature, contamination, and cross contamination of two-part materials. Recap immediately after use.

**Storage:** Keep containers tightly closed. Store product in original container at cool room temperature (< 25°C) and in a dry, well ventilated area. Avoid direct, strong light, sources of ignition and extremes of temperature. Shelf life for unopened product is two years from date of manufacture, provided that the material has been stored properly.

Specific uses: Dental material

#### SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION

**Eye protection, hands protection, Engineering controls or other equipment:** None required

#### SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES

**Color / Physical state:** Natural tooth shade or off-white thin paste.

**Odor:** Faint, characteristic

**Vapor pressure:**< 1 mm HgPH

**Solubility in water:** Nil

**Specific gravity:** 1.575

#### SECTION X - STABILITY AND REACTIVITY

**Conditions to avoid:** Heat (> 38°C), light, contamination.

**Materials to avoid:** Reducing and oxidizing agents, peroxides, amines.

**Hazardous decomposition products:** Under fire conditions and with amounts greater than that supplied in this product, hazardous polymerization may occur with heat build up and release of carbon monoxide, carbon dioxide, and oxides of nitrogen.

**Further information:** Stable material if stored and used as directed. Polymerization will occur when light-cured material is exposed to direct light.

#### SECTION XI - TOXICOLOGICAL INFORMATION

**Acute toxicity:** Not toxic. Minimal health hazard under normal conditions of use.

**Local effects:** Eye irritant, mucous membranes or skin on contact or with prolonged exposure. Prolonged or frequent skin contact may cause allergic skin reactions in some susceptible individuals.

#### SECTION XII - ECOLOGICAL INFORMATION

Avoid release of uncured material into the environment. To the best of our knowledge, polymerized material is inert. No other information is available. Follow all government regulations.

#### SECTION XIII - DISPOSAL CONSIDERATIONS

Dispose in accordance with federal, state and local regulations

#### SECTION XIV - TRANSPORT INFORMATION

DOT- Not regulated

#### SECTION XV - REGULATORY INFORMATION

Champ Pit and Fissure Sealant is a Class IIA Medical Device under the Medical Devices Directive 93/42/EEC.

<b>Hazard symbol:</b>	Xi
<b>Hazardous ingredient:</b>	Uncured acrylate ester monomers
<b>Risk phrases</b>	R36/37/38: Irritating to eyes, respiratory system and skin.
<b>R43:</b>	Sensitizing by skin contact
<b>Safety phrases:</b>	S24/25:Avoid contact with skin and eyes.
<b>S3/9/49:</b>	Keep only in the original container in a cool, well ventilated place.

#### SECTION XVI - OTHER INFORMATION

List of relevant R phrases: R36/37/38: Irritating to eyes, respiratory system and skin.
R43: Sensitizing by skin contact.

Training advice: Champ Pit and Fissure Sealant is for use by dental professionals only.

Restrictions on use: Champ Pit and Fissure Sealant is to be sold to and used by dental professionals only.
Further information: The information presented herein is believed to be factual as it has been derived from the works of persons believed to be qualified experts. However, nothing contained in this information is to be taken as a warranty or representation for which Centrix, Inc. bears legal responsibility. The user should review any recommendations in the specific context of the intended use to determine whether they are appropriate.
Sources of key data: National Institute for Occupational Safety (NIOSH) US Occupational Safety and Health Administration (OSHA)Eur-Lex European Union Law. Information which has been added, deleted or revised: This Safety Data Sheet has been revised to meet the requirements of Regulation (EC) No. 1907/2006 (REACH).

<b>centrix</b> <sup>®</sup> MAKING DENTISTRY EASIER	
 <b>CENTRIX INCORPORATED</b> <b>770 River Road, Shelton, CT 06484 USA</b> <b>(203) 929-5582 • (800) 235-5862 • FAX (203) 929-6804</b> <b>www.centrixdental.com</b>	
<b>EC REP</b> <b>Dental-Contact Vertriebs OHG, Im Spiet 101, DE 26506 Norden</b>	
<b>EMERGO AUSTRALIA</b> <b>Level 20, Tower II, Darling Park, 201 Sussex Street</b> <b>Sydney, NSW 2000, Australia</b>	

©2021 Centrix, Inc.	Date Revised: 07/2021	012086 REV F
---------------------	-----------------------	--------------