Encore® D/C & Encore® MiniMix® D/C Automix

ENGLISH

ENCORE® D/C CORE BUILDUP MATERIAL IN MINIMIX® & AUTOMIX DELIVERY SYSTEMS

DUAL-CURE COMPOSITE RESIN

Indicated for use as a core buildup composite under indirect restorations and as a post cement when used in conjunction with a suitable dentin adhesive.

INSTRUCTIONS FOR USE

TOOTH/CANAL PREPARATION

- Isolate the tooth to be restored. We recommend using a rubber dam or cotton rolls.
- Remove existing restorations and decay. We recommend using Expose® Caries Indicator to confirm that all carious material has been removed; if using a selective caries removal technique, use Expose to confirm that healthy dentin completely surrounds the carious lesion.
- If a post is to be used, prepare the post space according to the manufacturer’s directions.
- Acid etch the entire preparation for 20 seconds using 32% 40% phosphoric acid. We recommend using Centrix Onyx™ 40% phosphoric acid etching gel.
- Rinse with water spray for 30 seconds and lightly air dry. The enamel of the preparation should have a frosted appearance.
- Place the bonding system of choice into the preparation according to the manufacturer’s instructions. To ensure complete cure of all material in the canal, we recommend using Connexio®, a dual-cure 7th-generation adhesive that is a light-activated self-etch system with a dual-cure activator. Do not use self-etching, one-bottle bonding systems. When placing the primer, be sure that all of the solvent has been evaporated and the preparation has a glossy appearance. Before light curing, the bonding agent should be thinned with a light air spray.

CORE BUILDUP USING ENCORE D/C

- Select the desired shade of Encore D/C core buildup material
 - Natural Shade (A2) for use under all–ceramic restorations.
 - Contrast Shade (blue) to distinguish core material from tooth structure.
- Inject the material into the prep. Fill from the bottom up to reduce the chance of voids. Keep tip embedded in the material to avoid air entrapment.
 - When using Encore D/C MiniMix cartridges
 - Insert MiniMix cartridge into a Centrix Snap-Fit™ Syringe (fits any composite gun).
 - Remove cap from MiniMix cartridge and attach mixing nozzle.
 - Always purge a small amount from the tip to assure proper mixing and flow.
 - When using Encore D/C Automix 25 mL cartridges
 - Place the cartridge into a Type S25 Applicator Gun.
 - Remove the cap from the cartridge and attach a new mixing tip.
 - Always purge a small amount from the tip to assure proper mixing and flow. The Centrix Access Automix Needle Tip (REF 290025) may be snapped into the open end of the mixing tip to place material into the open canal.
- Cure material.
 - Light Cure – follow the recommendations for your light cure unit; units vary widely for both output and light spectrum used. Generally –
 - for a 2 mm depth of cure 20 seconds (:20”)
 - for a 4mm (max.) depth of cure 40 seconds (:40”)
 - Self Cure (Encore D/C is a dual cure material, it will completely self cure) –
 - Working Time – 1:30 minutes
 - Total Setting Time – 5:00 minutes
- Once set, cut and finish the preparation as normal using water-cooled diamonds or burs.
- To keep material fresh and avoid pre-mature setting, keep used mixing tip on cartridges between uses.

NOTE: When using a Matrix Band or prefabricated Core Form, always place the mixing tip directly into the preparation/matrix/core form and express the paste from the bottom up. Test the core material with an explorer to insure that material has set prior to removal of the matrix band. Hold the loaded core form onto the prepared tooth when light curing until the paste is fully set.

POST CEMENTATION USING ENCORE D/C

- Prepare post space per the manufacturer’s recommendations and as described above. We recommend using a fiber/ resin post to create a fully bonded ‘monoblock’ restoration that bonds together the tooth/post/core.
- Use a dual-cure bonding system (e.g., Connexio) to ensure that all material, even at the deepest part of the post space is fully cured. Follow the manufacturer’s instructions.

PRECAUTIONS

- If pulp is exposed, use a suitable liner. Do not use Encore D/C core buildup paste for direct application on dentin near the pulp or if the pulp chamber has been exposed.
- Contact with skin and eyes should be avoided. If accidental contact occurs, rinse immediately with plenty of water and consult a physician if necessary.
- Do not use with materials containing eugenol. This could inhibit polymerization.
- Water and oil-contaminated air must be avoided, as they inhibit polymerization at the contact area.
- Avoid contact with gingival tissue and adjacent teeth.
- Placement of any composite is contraindicated if a dry field is not possible.
- There are no known systemic side effects; however, do not use in patients with known sensitivities to Bis-GMA, benzoyl peroxide or composite resin materials.

STORE AWAY FROM SUNLIGHT, AT ROOM TEMPERATURE 23°C (72°F)

310150	Encore® DC, Natural, With Fluoride (25 mL)
310152	Encore® DC, Contrast Blue With Fluoride (25 mL)
310254	Encore® D/C MiniMix®, Natural, 36 (.5 mL)
310255	Encore® D/C MiniMix®, Contrast 36 (.5 mL)

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310130 Encore® D/C MiniMix®, Contrast 12 (.5 mL)

700000 Encore® D/C MiniMix®, Natural 12 (.5 mL) Kit:

Includes Snub Nose Snap-Fit Syringe

SAFETY DATA SHEET (SDS) | In accordance with Hazard Communication Standard (HCS), 29CFR 1910.1200(g), 2012 revision

SECTION I – PRODUCT/ COMPANY IDENTIFICATION

Product Trade Name: ENCORE DC MINIMIX

Part (Item) Number: See Above

Product Name: Dental Core Build-Up material

Manufacturer: CENTRIX, INC. Address: 770 River Road

City, State, Zip: Shelton, CT 06484 U.S.A.

Emergency or Information Telephone Numbers: (203)929-5582,

Toll-Free: (800)235-5862 (USA & Canada) Fax: (203)929-6804

Date Revised: 11/2020

SECTION II – HAZARDS IDENTIFICATION

Designation of Risk: Xi, irritant

Special Details of Risk (R-categories): 36-38

Detrimental effects/symptoms: Irritant in case of eye, skin and inhalation contact.

Harmful if swallowed. May cause sensitization by skin contact.

SECTION III – HAZARDOUS INGREDIENTS / IDENTITY INFORMATION

Hazardous Components:	%	CAS NO.	SYMBOL	R	OTHER LIMITS
Blend of multifunctional methacrylates	26 - 32	N/A	Xn/Xi	22/36/37/38	
Initiator	0.10 - 0.35	94-36-0	E/Xn	2/7/36/43	

Exact % content withheld due to trade secret restrictions.

SECTION IV – FIRST AID MEASURES

Skin: Immediately wash skin with soap and water. If irritation occurs seek

medical advice.

Ingestion: If swallowed do not induce vomiting. Give large quantities of water and milk. Seek

medical attention.

Inhalation: Remove to fresh air. If breathing becomes difficult, contact a physician.

Eyes: Flush with plenty of water. Seek immediate medical attention by an

Ophthalmologist.

SECTION V – FIRE-FIGHTING MEASURES

The product is not considered to be a fire hazard. The information below is given

for extinguishing a fire generated by other sources. Suitable extinguishing media: Water, Chemical

Foam, Carbon Dioxide, or Dry Chemical.

Forbidden extinguishing media: Not determined.

Special firefighting measures: None. Wear self-contained breathing apparatus.

Unusual fire and explosion hazards: Heat can cause polymerization with rapid release of energy.

Special protection equipment: Self Contained Breathing Apparatus if toxic fumes are generated.

SECTION VI – ACCIDENTAL RELEASE MEASURES

Personal precautions: Adopt some precautions of section VII.

Environmental precautions: Keep spilled material out of sewers.

Reclaiming methods: Absorb spills with inert material - keep spilled material out of sewers.

SECTION VII – Handling and Storage

Handling: Avoid unnecessary exposure. Storage: Room temperature. Containers:

Supplied by manufacturer.

Incompatibilities: Avoid the contact with reducing and oxidizing agents, peroxides and amines.

Other precautions: Use according to directions. Avoid contact with eyes, skin and soft tissue.

SECTION VIII – Exposure Controls / Personal Protection

General measures: Local Exhaust Ventilation: Good general ventilation should be sufficient to

control airborne levels of vapors released by uncured material.

Limits: TWA: Not determined; TLV: Not determined.

Respiratory protection: Avoid prolonged or excessive breathing of vapors of uncured material.

Eye protection: Safety glasses or goggles.

Skin protection: Respect good personal hygiene practice. Avoid unnecessary exposure.

Hands protection: Impervious gloves.

Other equipment: Lab coat.

SECTION IX – Physical and Chemical Properties

Appearance: Colored/uncolored paste Explosion limits: Not determined Odor: None

Water solubility: Insoluble Boiling point: Not determined Evaporation rate: Not determined

Flash point: Not determined Other Data: Not available or none.

SECTION X – Stability and Reactivity

Stability: Stable if stored as directed.

Conditions to avoid: Heat, light, aging and contamination.

Materials to avoid: Reducing and oxidizing agents, peroxides and amines.

Hazardous decomposition products: Oxides of Carbon.

Hazardous polymerization products: May occur.

SECTION XI – Toxicological Information

Ingestion: Uncured material may be harmful or fatal if swallowed.

Eye contact: May cause irritation and damage if not removed promptly.

Skin contact: Prolonged or repeated exposure to uncured material may cause irritation

or skin rash.

Inhalation: Prolonged or excessive inhalation may cause respiratory tract irritation.

Prolonged exposure: Not applicable.

SECTION XII – Ecological Information

There is no specific information on the product. Use according to good working practice.

Avoid releasing the product into environment.

SECTION XIII – Disposal Considerations

Unpolymerized (uncured) material may be hazardous waste.

Incinerate uncured material and dispose in accordance with local regulations.

SECTION XIV – Transport Information

Additional Information: Not dangerous cargo.

SECTION XV – Regulatory Information

CLASSIFICATION IN ACCORDANCE WITH EEC DIRECTIVES.

Hazard Classification Symbol: Xn

R22 Harmful if swallowed.

R36-38 Irritating to eyes, respiratory system and skin.

S36 Wear suitable protective clothing.

Ingredients determining the label: Blend of multifunctional methacrylate

SECTION XVI – Other Information

Product handling is exclusively by dental specialists per instructions for use.

CAUTION

The above information is based on presently available data and to our best knowledge

for handling the product under normal conditions. Any use of this product in any way not

indicated on this document or using it together with any other process/procedure will be

exclusively under the user’s responsibility.

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EMERGO AUSTRALIA

Level 20, Tower II, Darling Park, 201 Sussex Street

Sydney, NSW 2000, Australia