

REF 310108 72 Cubes

1

2

3-4

5

ENGLISH INSTRUCTIONS FOR USE

1. Peel backing off bottom of LubeCube.
2. Place and/or stick LubeCube in desired location (back of hand, bib, tray etc).
3. Insert tip of composite instrument into slit of LubeCube.
4. Glide tip of composite instrument back and forth in slit to clean and lubricate.
5. Remove composite instrument and begin sculpting the restoration.
6. Repeat steps 3 thru 5 as desired until restoration is complete.
7. Dispose of used LubeCube.

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SAFETY DATA SHEET (SDS)

In accordance with Hazard Communication Standard (HCS),
29CFR 1910.1200(g), 2012 revision

SECTION I - PRODUCT/ COMPANY IDENTIFICATION

PRODUCT NAME: LUBECUBE COMMON OR CHEMICAL NAME: Light-Cure Resin
Recommended Use: Dental Instrument Lubricant
PRODUCT (ITEM) NUMBER: 310108
MANUFACTURED BY: CENTRIX INC
EMERGENCY & INFORMATION TELEPHONE: (800) 235-5862 Address: 770 RIVER ROAD
City, State, Zip: SHELTON, CT 06484
GENERAL INFORMATION: 203-929-5582

SECTION II - HAZARDS IDENTIFICATION

For potential health hazards: SEE Section 11 below.
No other significant hazards to man or the environment are reasonably foreseen.

SECTION III - HAZARDOUS INGREDIENTS/IDENTITY INFORMATION

Product Contents: Light-Cure Bonding Agent Liquid
Hazardous Ingredients: None listed.
Exact % content withheld due to trade secret restrictions.

SECTION IV - FIRST AID MEASURES

Inhalation: No harm anticipated. Skin Contact: Wash with soap and water.
Eye Contact: Flush with copious amounts of water.
Ingestion: Induce vomiting if large amount of material

SECTION V - FIRE-FIGHTING MEASURES

Suitable Extinguishing Media: Water contact acceptable.
Extinguishing Media which should not be used: None.
Special hazards from combustion products: None specific.
Special protective equipment for fire-fighters: None.

SECTION VI - ACCIDENTAL RELEASE MEASURES

Personal Precautions: Avoid contact of skin or eyes.
Wash hands with soap and water if contact occurs.
Environmental Precautions: Dispose in accordance with all federal, state, and local regulations.
Methods for Cleaning Up: Wipe up with soap, water, and paper towels.

SECTION VII - HANDLING AND STORAGE

Handling Precautions: Not applicable.
Safe Storage Conditions: Store at temperatures not exceeding 75°F (24°C).
Other Restrictions or Recommendations: Avoid exposure to moisture and light.

SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION

Respiratory Protection & Ventilation: Not required.
Protective Clothing or Equipment: Not required under normal handling.
Unusual Work/Hygienic Practices: Avoid prolonged contact with skin; wash with soap and water.

SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Clear Liquid. Odor: Slight acrylic odor. pH: Not applicable.
Boiling Point/Range: Not applicable. Melting Point/Range: Not applicable.
Flash Point: Not applicable. Auto-Ignition Temperature: Not applicable. Flammability: No.
Explosive Limits in Air: Not applicable
Oxidizing Properties: None.
Vapor pressure (mm Hg): Negligible at room temperatures.
Vapor Density (Air=1): >1 Specific Gravity (H2O=1): ~1
Solubility in Water: Insoluble. Reactivity in Water: Not reactive.
Other Data: None or not available

SECTION X - STABILITY AND REACTIVITY

Stability: Generally stable.
Conditions of Instability: Avoid exposure to light.
Incompatibility with Various Substances: None
Hazardous Decomposition Products: None
Hazardous Polymerization: May occur upon exposure to light where large quantities of material are involved.
Conditions to Avoid: Avoid exposure to light and temperatures above 80°F.

SECTION XI - TOXICOLOGICAL INFORMATION

Potential Acute Health Effects: Possibility of skin irritation or sensitization if not washed off promptly.
Potential Chronic Health Effects: Possibility of skin sensitization or irritation upon prolonged or repeated exposure.
Signs and Symptoms of Exposure: Skin/eye irritation.
Medical Conditions Generally Aggravated by Exposure: None known.
Carcinogenic or Potential Carcinogenic Effects: No ingredient listed as a carcinogen.
Other Potential Health Effects: None known.

SECTION XII - ECOLOGICAL INFORMATION

No ecological damage or impact anticipated when handled as indicated.

SECTION XIII - DISPOSAL CONSIDERATIONS

Waste Disposal: Do not dispose into sink.
Dispose of in accordance with all federal, state, and local regulations.

SECTION XIV - TRANSPORT INFORMATION

UN Classification: Not applicable.
Identification: Not applicable.
Special Provisions for Transport: Avoid elevated temperatures and high intensity light.

SECTION XV - REGULATORY INFORMATION

U.S. FDA Dental Device: Class I product
EC Medical Device Directive: Class I product

SECTION XVI - OTHER INFORMATION

Product handling exclusively by dental specialists per instructions for use.

CAUTION

The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.

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EMERGO AUSTRALIA

Level 20, Tower II, Darling Park, 201 Sussex Street
Sydney, NSW 2000, Australia