# NoCord™ VPS

One Step, Self-Retracting, Impressioning System

## ENGLISH

### INSTRUCTIONS FOR USE — For Professional Use Only

The NoCord™ VPS Introductory Kit is a revolutionary one step, self-retracting, impressioning system designed specifically for the wash/tray technique. Trial consists of the Centrix NoCord™ Wash and NoCord™ MegaBody™ Tray Materials and all the accessories needed (trays, mixers and mixing tips) to provide indirect dental impressions. Both NoCord™ Wash and NoCord™ MegaBody Tray Materials are addition-reaction silicone-based elasticomeric materials, or vinyl polysiloxane (VPS), offering excellent physical and dimensional properties as well as providing high accuracy and stability. Both materials are color differentiated for function, with work and set times for each compatible with this technique.

Our unique NoCord™ Wash Material is so named because no retraction cord and/or astringent agent is needed to prepare the impression site in advance; NoCord™ Wash Material is a low viscosity, light body VPS material which retracts the gingiva and stops fluid flow when the compressive action of the high viscosity, MegaBody™ Tray Material is applied. For optimal margins without a separate retraction procedure, NoCord™ Wash and MegaBody Tray Materials have been specially formulated to work together as a system and should always be used together.

NoCord™ Wash Material also includes 15% alum (ammonium aluminum sulfate) as an astringent agent to control bleeding as a further measure to keep the site dry and to save steps and time. Having exposed the sub-gingival margin, NoCord™ Wash Material then cures with MegaBody™ Tray Material to provide a high detail, low shrinkage, clear, final impression.

**NOTE: For best results, always use NoCord VPS at room temperature. Use of NoCord VPS outside of recommended temperature range may affect set times.**

**INDICATIONS:** NoCord™ VPS system is designed to take impressions without separate gingival retraction and hemostatic control steps. Using the wash/tray technique, impressions can be achieved for temporary or permanent crowns or bridges, implants, veneers, aligners and denture cases.

**CAUTION:** NoCord™ Wash Material contains 15% alum (ammonium aluminum sulfate) as an astringent agent; patients sensitive to this astringent should not be exposed to this product. Astringents are irritants, flush with water immediately if tissue contact outside the oral cavity occurs.

#### CONTRAINDICATIONS:

- Do not use NoCord™ Wash Material with any heavy body tray material other than MegaBody™ Tray Material, or product performance may be affected.
- It is not necessary to use retraction cord or paste materials for most cases.
- For best impression results:
  - Dry prep with dressing or air syringe
  - A walled, dual arch or full arch tray is recommended such as Centrix 3.n.One™ Impression trays
- NoCord VPS materials (NoCord Wash & NoCord MegaBody Tray) should not come in contact with sources of sulfur (latex gloves, ferric sulfate solutions, etc.); setting may be affected. Store at room temperature, not to exceed 75°F (24°C).

**MODEL POURING:** For best results, NoCord VPS models should be poured 2-3 hours after removal from the mouth to allow the VPS material to degas. NoCord VPS remains dimensionally stable up to 20 days.

#### WORKING TIMES:

Product/Material	Color	Technique	Mixing Time	Working Time	Total Set Time
NoCord™ Wash	Coral Orange	Wash & Tray	Auto	1:00	4:45
NoCord™ MegaBody™ Tray	Seafoam Green	Wash & Tray	Auto	1:00	4:45

#### PROCEDURAL STEPS:

- Prepare tooth as usual.
- Dry prep with dressing or air syringe.
- Syringe NoCord Wash Material (Coral Orange) into the gingival sulcus and around the complete preparation. Use copious amount of the wash material.
- Fill the impression tray (a walled, dual arch or full arch tray) with NoCord MegaBody Tray Material (Seafoam Green).
- Place over the unset wash material and have patient bite down until the material is set.
- Remove from patient's mouth. Full set is 4:45 minutes.
- Rinse site with water tray removal to wash away any remaining excess material or astringent.
- Check with magnification to ensure all of the margins have been captured.

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