# Onyx™

## INSTRUCTIONS FOR USE

40% Phosphoric Acid Etching Gel

### GB - ENGLISH INSTRUCTIONS FOR USE

**Product Description:** Onyx™ is a viscous, high-contrast 40% phosphoric acid etchant (H3PO4). The black silicon carbide particles provide easier visibility for placing the etchant and then helps verify that the etchant has been cleanly and completely rinsed away.

- Prepare tooth using pumice, disc or diamond bur on uncut enamel. Rinse and dry thoroughly.
- Remove cap from syringe and screw on a disposable needle tip.
- Apply Onyx etching gel to enamel and/or dentin.
- Let stand on tooth surface for 15 seconds. Rinse with water and high-volume suction.
- Dry gently. Tooth surface should appear dull. If not, repeat above steps.
- Remove needle tip, re-cap and disinfect syringe barrel.

**INDICATIONS:** Appropriate for any bonding system requiring selective or total etching.

#### CAUTION & WARNINGS:

- Onyx contains 40% phosphoric acid; if in contact with any soft tissue, IMMEDIATELY rinse area with copious amounts of water.
- Verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check.
- Re-cap, disinfect and wipe syringe with an intermediate level disinfectant between uses.
- Dispose of tip and empty syringe properly.
- Near pulp exposures should be treated with a protective base prior to placing etchant.
- Do not use after expiration date noted on container.

**DISCLAIMER:** Centrix products have been developed for use in dentistry. We have provided written, verbal and online information as to the use of each product. This information does not relieve you from examining our products upon receipt to determine their suitability for use in the intended manner of application. The usage of our products is beyond our control and the user maintains full responsibility. We guarantee the quality of our products in accordance with the existing standards and corresponding to the conditions stated in our instructions for use and safety data sheets.

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### SAFETY DATA SHEET (SDS)

In accordance with Hazard Communication Standard (HCS), 29CFR 1910.1200(g), 2012 revision

**SECTION I - PRODUCT / COMPANY IDENTIFICATION**
Product Trade Name: Onyx, Etching Gel
Synonyms: 40% Phosphoric Acid etching gel
Part (Item) Number: See Above
Application: Etching gel for dental restorations
Manufacturer: Centrix Inc.
Address: 770 River Road
City, State, Zip: Shelton Ct, 06484
Emergency Telephone Number: 203-929-5582
Information Telephone Number: 800-235-5862

**SECTION II – HAZARD IDENTIFICATION**
**Classification of the Substance or Mixture:**
Acute toxicity, Oral (Category 4)
Acute toxicity, Dermal (Category 5)
Serious eye damage (Category 1)

**Risk Phrases:**
Symbol: C
**R34:** Causes burns.
**Label Elements:**
**Trade Name:** PHOSPHORIC ACID
**Hazard Statements:**
H302: Harmful if swallowed.
H313: May be harmful in contact with skin.
H314: Causes severe skin burns and eye damage.
H330: Fatal if inhaled.
**Precautionary Statements:**
P260: Do not breathe dust / fume / gas / mist / vapors / spray.
P280: Wear protective gloves / protective clothing / eye protection / face protection.
P284: Wear respiratory protection.
P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes.
Remove contact lenses if present and easy to do – continue rinsing.
P310: Immediately call a POISON CENTER or doctor / physician

**SECTION III - COMPOSITION / INFORMATION ON INGREDIENTS**
**Hazardous Components:** Phosphoric Acid Concentration: 40%
**CAS No.:** 7664-38-2
**EC Number:** 231-633-2
**Index Number:** 015-011-00-6
**Molecular Weight:** 98 g/mol

**SECTION IV - FIRST AID MEASURES**

Skin: Immediately wash skin with plenty of water.
If irritation occurs seek medical advice.
Ingestion: If swallowed do not induce vomiting. Give large quantities of water and milk.
Seek medical attention.
Eyes: Flush with plenty of water.
Seek immediate medical attention by an Ophthalmologist.

**SECTION V - FIRE-FIGHTING MEASURES**

The product is not considered to be a fire hazard. The information below are given for extinguishing a fire generated by other sources.
Extinguish Extinguishing Media: Water, Chemical Foam, Carbon Dioxide, or Dry Chemical
Forbidden Extinguishing Media: None.
Special Exposure Hazards: Toxic fumes of Phosphoric Acid.
Unusual fire and explosion hazards: Contact with some metals may cause formation of hydrogen gas.
Special Protection Equipment: Self Contained Breathing Apparatus if toxic fumes are generated.

**SECTION VI - ACCIDENTAL RELEASE MEASURES**

Personal Precautions: See sections VII and VIII.
Environmental Precautions: None.
Methods for Cleaning up: Clean up using impervious gloves and dispose of in an approved manner.

**SECTION VII - HANDLING AND STORAGE**
Handling: Avoid unnecessary exposure.
Storage: Room Temperature.
Containers: Supplied by manufacturer.
Incompatibility: Strong caustics, sulfites, sulfides and metals.
Other Precautions: Use according to directions.
Avoid contact with eyes, skin and soft tissue.

**SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION**
General Measures Ventilation: N/A
Limits: TWA: 1mg/m3; TLV: 3 mg/m3.
Respiratory Protection: N/A
Eye Protection: Safety Glasses or Goggles.
Skin Protection: Respect good personal hygiene practice.
Avoid unnecessary exposure.
Hands Protection: Impervious gloves
Other Equipment: Lab Coat.

**SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES**
Appearance: Gel Explosion Limits: N/A
Odor: None
pH 10% solution: below 2
Boiling Point: 100°C
Flash Point: N/A

Water solubility: Soluble.

Specific Gravity: 1.3 g/cm3

Evaporation Rate: N/A

**SECTION X - STABILITY AND REACTIVITY.**

Stability: Stable.
Conditions to avoid: Temperatures over 40°C
Materials to avoid: Strong Caustics, Sulfites, Sulfides, Metals.
Hazardous Decomposition Products: Toxic fumes if temperature is above 100°C.
Hazardous Polymerization Products: None

**SECTION XI - TOXICOLOGICAL INFORMATION.**
Ingestion: May cause sore throat, abdominal pain, nausea and severe burns of mouth, throat and stomach.
Eye contact: May cause redness, pain, blurred vision, eye burns or permanent eye damage.
Skin Contact: May cause redness, pain and severe skin burns.
Inhalation: Inhalation is not a hazard unless the product is heated at high temperatures.
In these conditions may cause irritation of the respiratory system, asphyxia or permanent damage to the lungs.
Prolonged Exposure: Burns to eye, soft tissue or respiratory tract.

**SECTION XII - ECOLOGICAL INFORMATION.**
There is no specific information on the product. Use according to good working practice. Avoid releasing the product into environment.

**SECTION XIII - DISPOSAL CONSIDERATIONS.**
Dispose in accordance with local authorities.

**SECTION XIV - TRANSPORT INFORMATION.**
Substances determining Hazard: Phosphoric Acid.
SEA TRANSPORTATION: IMCO # 8204 UN # 1805
AIR TRANSPORTATION: ICAO/IATA CLASS: 8
GROUND TRANSPORTATION: RID/ADR CLASS: 9
SUBSTANCE IDENTIFICATION: # 1805 ITEM # 11C
HAZARD LABEL: 9
HAZARD IDENTIFICATION: 80

**SECTION XV - REGULATORY INFORMATION**
Hazard Classification Symbol: C
R34 Causes burns
S24/S25 Avoid contact with skin and eyes
S26/S28 In case of contact with eyes or skin, rinse immediately with plenty of water and seek medical advice.

**SECTION XVI - OTHER INFORMATION**
Product handling exclusively by dental specialists per instructions for use.

**CAUTION:**
The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.

	<b>CENTRIX INCORPORATED</b> <b>770 River Road, Shelton, CT 06484 USA</b> <b>(203) 929-5582 • (800) 235-5862 • FAX (203) 929-6804</b> <b>www.centrixdental.com</b>
	 <b>Dental-Contact Vertriebs OHG, Im Spiet 101, DE 26506 Norden</b>

<b>EMERGO AUSTRALIA</b> <b>Level 20, Tower II, Darling Park, 201 Sussex Street</b> <b>Sydney, NSW 2000, Australia</b>		
<b>©2021 Centrix, Inc.</b>	<b>Date Revised: 10/2021</b>	<b>000835 REV K</b>