REF 390010

REF 390011

ENGLISH • INSTRUCTIONS FOR USE

ProxiDiscs[®] For Proximal Separation

REF 390010 – ProxiDiscs Diamond (8 µm), 1-sided, 4 pk
REF 390011 – ProxiDiscs Smooth, 2-sided, 8 pk

PRODUCT DESCRIPTION:

ProxiDiscs are stainless steel proximal separator discs designed to help adjust, break or smooth contacts, or to remove excess bonding agents, resins, or cements from interproximal areas. There are two ProxiDisc versions: Smooth (REF # 390011) enables easy separation of tight contacts; the diamond-coated ProxiDisc (REF # 390010) has a fine diamond coating on one side to help in adjusting tight contacts or removing excess materials from one tooth only, while the smooth surface on the opposite side protects the adjacent tooth.

INSTRUCTIONS FOR USE:

- 1) Choose the preferred disc, (smooth sides, or diamond coated where an abrasive effect is preferred).
- 2) Tie a safety loop of dental floss through the hole in the handle.
- 3) Using a firm rocking motion, roll the ProxiDisc through the contact area with pressure to remove excess material trapped in tight interproximal area.
- 4) If additional reduction is required, move the ProxiDisc buccolingually until you feel the restoration settling onto the adjacent tooth as the abrasive diamond coating relieves the pressure.
- 5) Clean and sterilize after each use.

MAINTENANCE AND STERILIZATION:

- ProxiDiscs are supplied non-sterile. Prior to first use, and immediately after each use, the instruments must be disinfected or sterilized.
 - Keep ProxiDiscs separate from non-stainless steel instruments, such as rubber polishers and abrasive stones.
 - Use only cleaning/disinfecting solutions that provide corrosion protection and strictly observe the concentrations recommended by the manufacturer. Rinse thoroughly with water after immersion into cleaning/disinfecting solution, otherwise discoloration or staining may occur.
- Sterilization:** For sterilization processes, we refer to the ISO 17664 standard and recommend the following procedures:
- | | |
|---|---|
| 1) 273°F/134°C Cycles: | 2) 250°F/121°C Cycles: |
| • T _{min.} = 134°C T _{max.} = 138°C | • T _{min.} = 121°C T _{max.} = 125°C |
| • Pressure = 45 PSI/3.15 bar (abs) | • Pressure = 30.5 PSI/2.10 bar (abs) |
| • Time = 4 minutes | • Time = 16 minutes |

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