^{CENTRIX} RapidCore[®] MiniMix[®]

Self-Cure Composite Resin Core Buildup Material

ENGLISH

RapidCore[™] MiniMix[™]

Self-Cure Composite Resin Core Buildup Material

REF 310126 Natural, 12 Units
REF 310256 Natural, 36 Units
REF 310127 White, 12 Units
REF 310257 White, 36 Units
REF 310128 Contrast, 12 Units
REF 310258 Contrast, 36 Units

Storage: Refrigeration may be used to extend shelf life. Return to room temperature before use.

Indications for use:

For dentin replacement and buildup under prefabricated crowns and bridges

1. Prepare the tooth:

- Remove all decay and old restorative material.
- Any eugenol present must be rinsed away with ethanol before proceeding.

2. Etch the tooth:

- Attach needle to syringe and express with a slow, steady pressure.
- Apply Onyx[™] 40% phosphoric acid etching gel to enamel first, then dentin.
- Wait 10 seconds, then rinse with water and dry gently.
- Surface should be moist, yet non-shiny.

3. Apply bonding agent:

- Any agent compatible with self-cure composites.

4. Build up core:

- Insert RapidCore Mini-Mix cartridge into any C-R[®] Syringe and remove cap.
- Attach mixing nozzle to the cartridge.
- Purge a small amount from the tip to assure proper mixing and flow.
- Inject the material into the prep, keeping the needle embedded in the material to reduce the chance of voids. fill from the bottom up. Working time is approximately 60 seconds.

5. Finish the core:

- After material is fully set (typically 2 minutes after injection), RapidCore can be cut to the desired shape, using water-cooled diamonds or burs.

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SAFETY DATA SHEET (SDS) In accordance with Hazard Communication Standard (HCS), 29CFR 1910.1200(g), 2012 revision

SECTION I - PRODUCT/ COMPANY IDENTIFICATION

Product Trade Name: RapidCore™ Self-Cure Composite Resin Core Buildup Material
Synonyms: Self-Cured Core Build-up Restorative
Recommended Use: Dental Material
Part (Item) Number:
REF 310126 Natural, 12 Units REF 310256 Natural, 36 Units
REF 310127 White, 12 Units REF 310257 White, 36 Units
REF 310128 Contrast, 12 Units REF 310258 Contrast, 36 Units

Manufacturer: CENTRIX INC
Address: 770 RIVER ROAD
City, State, Zip: SHELTON, CT 06484
Information Telephone Number: (800) 235 - 5862 (203) 929-5582
Emergency Telephone Number: (800) 235 - 5862

SECTION II – HAZARDS IDENTIFICATION

For potential health hazards: SEE Section 11 below.
No other significant hazards to man or the environment are reasonably foreseen.

SECTION III - COMPOSITION / INFORMATION ON INGREDIENTS

Product Contents: Self-Cured Core Build-up Restorative Material
Hazardous Ingredients: None listed.
Exact % content withheld due to trade secret restrictions.

SECTION IV - FIRST AID MEASURES

Skin: Wash with soap and water.
Ingestion: No harm perceived, however, if large quantities are accidentally swallowed, seek medical attention.
Inhalation: No harm anticipated from materials in quantities as packaged.
Eyes: Flush out with copious amounts of water. Seek medical attention if irritation persists.

SECTION V - FIRE-FIGHTING MEASURES

Suitable extinguishing media: Water contact acceptable.
Forbidden extinguishing media: None.
Special exposure hazards from combustion products: None in quantities as packaged.
Special protection equipment: None.

SECTION VI - ACCIDENTAL RELEASE MEASURES

Personal precautions: Avoid contact of skin or eyes.
Wash hands with soap and water if contact occurs.
Environmental precautions: Dispose in accordance with all local, state, and federal regulations.
Reclaiming methods: Wipe up and clean with soap and water.

SECTION VII - HANDLING AND STORAGE

Handling: Avoid exposure to moisture and light.
Storage: Store at temperatures not exceeding 73° F (23°C).
Other precautions: Avoid exposure to moisture and light.

SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION

General measures: Avoid prolonged contact with skin; wash with soap and water.
Respiratory protection: Not required.
Skin protection: Not required.
Hands protection: Not required.
Other equipment: Not required.

SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Part A: White, Natural or Blue Part B: Off-white paste
pH: All Components: Not Applicable. Flammability: No
Odor: Negligible Vapor pressure (mm Hg): Negligible at ambient temperatures.
Vapor Density (Air=1): Not Applicable. Boiling point: Not Applicable.
Specific Gravity (H2O=1): >1 Flash point: Not Applicable.
Auto-Ignition Temperature: Not Applicable. Melting Point/Range: Not Applicable.
Reactivity in Water: None. Explosive Limits in Air: Not Applicable.
Oxidizing Properties: None. Solubility in Water: Virtually insoluble.
Other Data: Not available or none.

SECTION X - STABILITY AND REACTIVITY

Stability: Generally stable unless exposed to elevated temperatures and high intensity light.
Conditions of Instability: Avoid exposure to elevated temperatures and high intensity light.
Incompatibility with Various Substances: None.
Hazardous Decomposition Products: None.
Hazardous Polymerization: Will not occur.
Conditions to Avoid: Elevated temperatures and high intensity light.

SECTION XI - TOXICOLOGICAL INFORMATION

Potential Acute Health Effects: Possibility of skin irritation or sensitization if not washed off promptly.
Potential Chronic Health Effects: Possibility of skin sensitization or irritation upon prolonged or repeated exposure.
Signs and Symptoms of Exposure: Skin/eye irritation.
Medical Conditions Generally Aggravated by Exposure: Previous skin sensitization.
Carcinogenic or Potential Carcinogenic Effects: No ingredient listed as a carcinogen.
Other Potential Health Effects: None known.

SECTION XII - ECOLOGICAL INFORMATION

No ecological damage or impact anticipated upon exposure to the environment in quantities as packaged.

SECTION XIII - DISPOSAL CONSIDERATIONS

Waste Disposal: Do not dispose into sink.
Dispose in accordance with all local, state, and federal regulations.

SECTION XIV - TRANSPORT INFORMATION

UN Classification: Not Applicable.
Identification: Not Applicable.
Special Provisions for Transport: Avoid elevated temperatures and high intensity light.

SECTION XV - REGULATORY INFORMATION


U.S. FDA Dental Device: Class II product.
EC Medical Device Directive: Class IIa product.

SECTION XVI - OTHER INFORMATION

Product handling exclusively by dental specialists per instructions for use.

CAUTION

The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.

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EMERGO AUSTRALIA

Level 20, Tower II, Darling Park, 201 Sussex Street
Sydney, NSW 2000, Australia

centrix®

Making Dentistry Easier.™