ENGLISH

33052X / 6 x 2g + 30 Needles; 33053X / 1 x 2g + 5 Needles

Light-cure flowable composite

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SAFETY DATA SHEET (SDS)

In accordance with Hazard Communication Standard (HCS), 29CFR 1910.1200(g), 2012 revision

SECTION I - PRODUCT / COMPANY IDENTIFICATION

Product Name: VERSAFLO®
Synonyms: Light-Activated Flowable Composite Resin
Recommended Use: Dental Material
Part (Item) Number: 330520, 330521, 330522, 330523, 330528,
330530, 330531, 330532, 330533, 330538
Company Name: Centrix Inc.
Address: 770 River Road
City, State, Zip: Shelton Ct, 06484
Emergency Telephone Number: 203-929-5582
Information Telephone Number: 800-235-5862

SECTION II - HAZARDS IDENTIFICATION

For potential health hazards: SEE Section XI below. 
No other significant hazards to man or the environment are reasonably foreseen.

SECTION III - COMPOSITION / INFORMATION ON INGREDIENTS

Product Contents: Flowable Light-Cure Composite Restorative / Veneer Cement
Hazardous Ingredients: None listed.
Exact % content withheld due to trade secret restrictions.

SECTION IV - FIRST AID MEASURES

Inhalation: No harm anticipated.
Skin Contact: Wash with soap and water.
Eye Contact: Flush immediately with copious amounts of water.
Ingestion: Induce vomiting if large amount of material is swallowed.

SECTION V - FIRE-FIGHTING MEASURES

Suitable Extinguishing Media: Water contact acceptable.
Extinguishing Media which should not be used: None.
Special exposure hazards from combustion products: None specific.
Special protective equipment for fire-fighters: None.

SECTION VI - ACCIDENTAL RELEASE MEASURES

Personal Precautions: Avoid contact of skin or eyes.
Wash hands with soap and water if contact occurs.
Environmental Precautions: Dispose of in accordance to all federal, state,
local regulations.
Methods for Cleaning Up: Wipe up with soap, water and paper towels.

SECTION VII - HANDLING AND STORAGE

Handling Precautions: None.
Safe Storage Conditions: Store at temperatures not exceeding 75°F (24°C).
Other Restrictions or Recommendations: Avoid exposure to moisture and light.

SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION

Respiratory Protection & Ventilation: Not required.
Protective Clothing or Equipment: Not required under normal handling.
Unusual Work/Hygienic Practices: Avoid prolonged contact skin;
wash with soap and water.

SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Thin paste. Odor: Negligible odor. pH: Not Applicable.
Boiling Point / Range: Not Applicable. Melting Point / Range: Not Applicable.
Flash Point: Not Applicable. Auto-Ignition Temperature: Not Applicable.
Flammability: No. Explosive Limits in Air: Not Applicable.
Oxidizing Properties: None. Vapor pressure (mm Hg): Negligible
Other Data: Not available or none.

SECTION X - STABILITY AND REACTIVITY.

Stability: Generally stable.
Conditions of Instability: Avoid exposure to light.
Incompatibility with Various Substances: None.
Hazardous Decomposition Products: None.
Hazardous Polymerization: May occur upon exposure to elevated
temperatures when large quantities of material are involved.
Conditions to Avoid: Avoid exposure to light and temperatures above 80°F.

SECTION XI - TOXICOLOGICAL INFORMATION.

Potential Acute Health Effects: Possibility of skin irritation or sensitization if not
washed off promptly.
Potential Chronic Health Effects: Possibility of skin sensitization or
irritation upon prolonged or repeated exposure.
Signs and Symptoms of Exposure: Skin/eye irritation.
Medical Conditions Generally Aggravated by Exposure: None known.
Carcinogenic or Potential Carcinogenic Effects: No ingredient listed as a
carcinogen.
Other Potential Health Effects: None known.

SECTION XII - ECOLOGICAL INFORMATION.

No ecological damage or impact anticipated when handled as indicated.

SECTION XIII - DISPOSAL CONSIDERATIONS.

Waste Disposal: Do not dispose into sink.
Dispose of in accordance with all federal, state, and local regulations.

SECTION XIV - TRANSPORT INFORMATION.

UN Classification: Not Applicable. Identification: Not Applicable.
Special Provisions for Transport: Avoid elevated temperatures and
high intensity light.

SECTION XV - REGULATORY INFORMATION

U.S. FDA Dental Device: Class II product
EC Medical Device Directive: Class IIa product

SECTION XVI - OTHER INFORMATION

Product handling exclusively by dental specialists per instructions for use.

CAUTION: The above information is based on presently available data and to our
best knowledge for handling the product under normal conditions. Any use of this
product in any way not indicated on this document or using it together with any
other process/procedure will be exclusively under the user's responsibility.

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EMERGO AUSTRALIA

Level 20, Tower II, Darling Park, 201 Sussex Street
Sydney, NSW 2000, Australia