### English

3306xx 20 x 0,25g

Light-cure hybrid composite restorative. Use with a light-cure bonding agent.

Do not use in presence of eugenol!

- 1) Prepare the tooth. Etch the surface, and apply a light-cure bonding agent (Adhere™).
- 2) Select a tip with the proper VITA® shade. Remove the sealing cap, and insert into a CENTRIX C-R® Syringe handle (NuGUN®, S/Class®, or Mk IIIp™ recommended).
- 3) Squeeze syringe with a slow, steady pressure - do not force. Test on pad before using in the mouth.
- 4) Place the tip at the base of the restoration and inject, withdrawing slightly as you squeeze. Place in layers no greater than 2mm thick. Do not remove the air-inhibited layer between applications. Cure each layer with a minimum 30 second exposure from a curing light in good working order.
- 5) Shape and finish, using standard composite finishing instruments. Polish using CENTRIX Min-Identoflex® polishing cups and points.

Remove VERSALITE™ tip from syringe and dispose of properly after use. Do not re-use the tip on another patient. Sterilize syringe gun by appropriate method. Store VERSALITE™ in a cool, dry place, not to exceed 75°F (24°C). Keep away from sunlight and sources of eugenol. Certain individuals may be sensitive or allergic to composite resins.

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# SAFETY DATA SHEET (SDS)

In accordance with Hazard Communication Standard (HCS), 29CFR 1910.1200(g), 2012 revision

## SECTION I - PRODUCT / COMPANY IDENTIFICATION

Product Name: VERSALITE®  
Synonyms: Light-Activated Hybrid Composite Resin  
Recommended Use: Dental Material  
Part (Item) Number: 330600, 330601, 330602, 330604  
330605, 330606, 330609, 330610  
Company Name: Centrix Inc.  
Address: 770 River Road  
City, State, Zip: Shelton CT, 06484  
Emergency Telephone Number: 203-929-5582  
Information Telephone Number: 800-235-5862

## SECTION II - HAZARDS IDENTIFICATION

For potential health hazards: SEE Section 11 below.  
No other significant critical hazards to man or the environment are reasonably foreseen.

## SECTION III - COMPOSITION / INFORMATION ON INGREDIENTS

Product Contents: Light-Cure Composite Restorative  
Components Containing Hazardous Ingredients: None listed  
Exact % content withheld due to trade secret restrictions.

## SECTION IV - FIRST AID MEASURES

Inhalation: No harm anticipated from materials in quantities as packaged.  
Skin Contact: Wash with soap and water.  
Ingestion: Induce vomiting if large amount of material is swallowed.  
Eyes Contact: Flush with copious amounts of water.

## SECTION V - FIRE-FIGHTING MEASURES

Suitable Extinguishing Media: Water contact acceptable.  
Forbidden Extinguishing Media which: None.  
Unusual Fire and Explosion Hazards: None in quantities as packaged.  
Special Protective Equipment: None.

## SECTION VI - ACCIDENTAL RELEASE MEASURES

Personal Precautions: Avoid contact of skin or eyes. Wash hands with soap and water.  
Environmental Precautions: Dispose of in accordance to all federal, state, local regulations.  
Methods for Cleaning Up: Wipe up with soap, water and paper towels.

## SECTION VII - HANDLING AND STORAGE

Handling Precautions: Not applicable.  
Safe Storage Conditions: Store at temperatures not exceeding 75°F (24°C).  
Other Precautions: Avoid exposure to moisture and light.

## SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION

Respiratory Protection & Ventilation: Not required.  
Protective Clothing or Equipment: Not required under normal handling.  
Unusual Work/Hygienic Practices: Avoid prolonged contact skin; wash with soap and water.

## SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Toothshade paste Odor: Negligible odor.  
pH: N/A Boiling Point/Range: N/A Melting Point/Range: N/A  
Flash Point: N/A Auto-Ignition Temperature: N/A Flammability: No  
Explosive Limits in Air: N/A Oxidizing Properties: N/A  
Vapor Pressure (mm Hg): Negligible at room temperatures.  
Vapor Density (Air=1): N/A Specific Gravity (H2O =1):~2.2  
Solubility in Water: Virtually insoluble.  
Reactivity in Water: Not reactive.  
Other Data: Not available or none.

## SECTION X - STABILITY AND REACTIVITY.

Stability: Generally stable.  
Conditions of Instability: Avoid exposure to elevated temperatures and light.  
Incompatibility with Various Substances: Not applicable.  
Hazardous Decomposition Products: Not applicable.  
Hazardous Polymerization: None in quantities as packaged.  
Conditions to Avoid: Exposure to light and temperatures above 75°F (24°C).

## SECTION XI - TOXICOLOGICAL INFORMATION.

Potential Acute Health Effects: Possibility of skin irritation or sensitization if not washed off promptly  
Potential Chronic Health Effects: Possibility of skin sensitization or irritation upon prolonged or repeated exposure.  
Signs and Symptoms of Exposure: Skin/eye irritation.  
Medical Conditions Generally Aggravated by Exposure: None known.  
Carcinogenic or Potential Carcinogenic Effects: No ingredient listed as a carcinogen.  
Other Potential Health Effects: None known.  
Prolonged Exposure: Burns to eye, soft tissue or respiratory tract.

## SECTION XII - ECOLOGICAL INFORMATION.

No ecological damage or impact anticipated upon exposure to environment in quantities as packaged

## SECTION XIII - DISPOSAL CONSIDERATIONS.

Waste Disposal: Do not dispose into sink. Dispose in accordance with all federal, state, and local regulations.

## SECTION XIV - TRANSPORT INFORMATION.

UN Classification: Not applicable.  
Identification: Not applicable.  
Special Provisions for Transport: Avoid elevated temperatures and high intensity light.

## SECTION XV - REGULATORY INFORMATION

U.S. FDA Dental Device: Class II product  
EC Medical Device Directive: Class IIa product

## SECTION XVI - OTHER INFORMATION

Product handling exclusively by dental specialists per instructions for use.

## CAUTION

The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.

Product covered by one or more of the following patents: US4,963,093 and US D353,673. Other U.S. and International patents pending or granted.

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### EMERGO AUSTRALIA

Level 20, Tower II, Darling Park, 201 Sussex Street  
Sydney, NSW 2000, Australia