

SuperCure and SuperCure-Q®

Light-Activated, Deep Curing Core Buildup Material

Core Buildup in 60 Seconds or Less with SuperCure®.

The easiest way to speed up your core buildups



SuperCure® Natural .25g Tip

SuperCure® Natural in EZ™ Gun

Single-dose. Great price. Problems solved.

With SuperCure and SuperCure-Q, no mixing or automix nozzles are required. Simply inject directly into the prep and build a high strength core in 60 seconds or less. SuperCure and SuperCure-Q are light-activated core buildup materials packaged in single-dose pre-filled tips. The non-slumping formula eliminates the need to use a matrix and the translucent material cures deep with any light source.

Features:

- Fast light-activation and high depth of cure - 8.9mm in 40 seconds
- Single-dose, prefilled syringe tips for maximum asepsis
- Non-slumping formula so there is no need to use a matrix
- No automix nozzles
- Available in two shades - natural and contrast
- Radiopaque

SuperCure – The original barium borosilicate glass filler that cuts like dentin.

SuperCure-Q – Reinforced formula of quartz class filler that cuts harder than dentin.



centrix®

Making Dentistry Easier.™

SuperCure and SuperCure-Q®

Light-Activated, Deep Curing Core Buildup Material

6-Step Guide To Building A Void-Free, Bonded, Composite Post & Core



1. Prepare the tooth. Bevel the edges to eliminate sharp angles.



2. Etch the post hole, rinse thoroughly and dry.



3. Apply bonding agent according to manufacturer's instructions.



4. Insert prefilled SuperCure tip to the bottom of the prep and inject a 2mm layer of material. Cure with light for 20 seconds.



5. Continue building up the core by injecting and slowly withdrawing the tip.



6. Light activate for 40 seconds and then cut and shape the core. The result is a solid, securely bonded post and core.

Original SuperCure and SuperCure-Q Specifications

	Original SuperCure	SuperCure-Q
Comprehensive Strength	235 Mpa	235 Mpa
Flexural Strength	87 Mpa	104 Mpa
Diametral Tensile Strength	50 Mpa	50 Mpa
Filler type	Barium borosilicate glass	Barium borosilicate and quartz glass
Cure depth, (ADA), 40 seconds	8.9mm	8.9mm



Natural



Contrast Blue



ORIGINAL SUPERCURE

Original formula of barium borosilicate glass filler cuts like dentin.

- REF 310140 30 Natural, .25g tip,
- REF 310142 30 Contrast Blue, .25g tip
- REF 310141 30 Natural, .50g tip
- REF 310143 30 Contrast Blue, .50g tip

Suggested Accessory

- REF 130650 Snap-Fit Gun



SUPERCURE-Q

Reinforced formula with quartz glass filler, cuts harder than dentin.

- REF 310146 30 Natural, .50g tip
- REF 310147 30 Contrast Blue, .50g tip

Suggested Accessory

- REF 130650 Snap-Fit Gun

Ordering is easy. Call 800.235.5862 or contact your preferred dealer. Visit centrixdental.com for more information.