



PARTNER



GROW



LEARN



CELEBRATING OVER 50 YEARS IN BUSINESS!

PARTNER, GROW & LEARN WITH US AS WE CONTINUE MAKING DENTISTRY EASIER FOR THE NEXT 50 YEARS.

As a leader in materials and delivery systems for over 50 years, Centrix manufactures and engineers hundreds of dental products focused primarily on general restorative and preventive dentistry, the heart and soul of most dental practices.

“Like most practicing dentists, I experienced problems with products that made the delivery of quality dentistry more difficult than it needed to be. My first products, AccuDose® tubes and the C-R® syringes, changed the use of dental materials and made dentistry easier. Now, over 50 years later, I am proud that Centrix continues to introduce products and ideas that fulfill our mission.” **Making Dentistry Easier. ~ Dr. William Dragan**



1970

1970

1984

1987

1989



TUBES AND PLUGS



MARK IIIP™ GUN & SNAP-FIT™ GUN



TEMPIT®



NEEDLETUBE®



BENDA® BRUSH

**FOR OUR FULL
INNOVATION
TIMELINE, VISIT:
[CENTRIXDENTAL.COM/
INDUSTRY-FIRSTS](http://CENTRIXDENTAL.COM/INDUSTRY-FIRSTS)**



PREVENTION FOR LIFE®

Caries rates in kids are decreasing! However, caries rates in adults are increasing, rapidly, as patients decline to pay for treatment not covered by insurance – or are not offered preventive care, as a matter of course.

The Centrix Prevention for Life program provides tools and information to help offices update their preventive treatment protocols and how they discuss them with patients.

Webinars and in-office Lunch & Learn programs that provide CE credits.

Risk Assessment Tools that help you quickly, easily and visually share with your patients your diagnosis. **Prescription and Treatment Program** that increases patient health (as well as practice health).

Products that are easy to place and your patients will love!

At Centrix, we believe that structured Risk Assessment quickly results in diagnosis of risk and a prescription for that patient (rinse, varnish, etc.), treatment plan and frequency for each patient, and easy patient acceptance!

START A PREVENTION FOR LIFE® PROGRAM TODAY!
VISIT CENTRIXDENTAL.COM/PFL
OR CALL 800.235.5862.



ADA Center for Evidence-Based Dentistry (November 2013) Topical Fluoride for Caries Prevention. Retrieved from <http://www.ada.org>

ORAL CARE AND PREVENTION – CENTRIX MAKES PREVENTION EASIER FOR YOU AND YOUR PATIENTS



LolloCaine®

Topical Anesthetic | 20% Benzocaine Gel

- 4 Great flavors: Bubble Gum, Cherry, Mint, Piña Colada.
- Rapid onset with no alcohol aftertaste!
- Unit-dose LolloPack® saves time and reduces cross-contamination.
- The LolloPack is the perfect “take-home pack” for post-op discomfort.



Exposé®

Caries Indicator on a Stick

- Easily identify and remove all infected dentin.
- Single-patient applicators make tray set-up easy.
- No syringe, no sheath, no tips to replace.
- No cross-contamination, no messy spills, no stains!
- Made in the USA.



FluoroDose®

Award-Winning Fluoride Varnish with Xylitol

- Tastes Better! Feels Better! Better Priced!
- Application takes less than a minute.
- Unit-dose delivery means easier use, easier cleanup.
- Available in 5 great-tasting flavors!



Caramel



Bubble Gum



Mint



Cherry



Melon



Champ®

Pit & Fissure Sealant | Glass-Filled Acrylate Resin

- **Bioactive:** hydrophilic resin forms a strong resin-hydroxyapatite complex that seals out micro-leakage.
- **Moisture-friendly formula:** tolerates moisture making it easier to bond to the tooth.
- **Margin Integrity:** lightly filled (37%) resin is durable and long-lasting; adapts to and seals the margins, helping to eliminate micro-leakage.

NOW AVAILABLE!

SilverSense SDF™

The 38% Silver Diamine Fluoride Sensation will Transform the Way You Treat Your Patients!

- Alleviates hypersensitivity
- Non-invasive procedure
- Application is fast, and easy
- Hardens infected dentin
- Cost-effective solution



LEARN MORE BY VISITING
CENTRIXDENTAL.COM/SILVERSENSE-SDF

TEMPORIZATION – COMPLETE SOLUTIONS DESIGNED TO WORK BETTER TOGETHER



Tempit® Moisture-Activated Temporary Filling & Sealing Material

- Ideal for endo access openings and short-term temporaries.
- Moisture-activated, sets in 5 minutes or less.
- Expands slightly when set, seals out microleakage and bacterial ingress.
- Contains calcium sulfate and zinc oxide.
- Significant bacteriocidal activity up to 14 days.¹

¹Antibacterial Properties of Temporary Filling Materials, JOE, Volume 32, Number 3, March 2006



Tempit®-E Moisture-Activated Temporary Filling & Sealing Material with Eugenol

- Palliative eugenol formulation for temporary filling use, especially in cases of reversible pulpal inflammation.
- Anti-inflammatory and local anesthetic effects helps facilitate pulpal healing.
- Moisture activated, sets in five minutes.
- Expands slightly when set to seal out microleakage and bacterial ingress.



Tempit® L/C Light-Activated Temporary Material for Inlays & Onlays

- Perfect for indirect restorations, inlays/onlays.
- Seals implant screw access openings.
- Sets firm yet flexible.
- Easily removed with an explorer without damaging the margins.
- Microfilled polymer resin.
- Light-activated, sets in 40 seconds or less.
- No etching or bonding agent required.



Tempit® Ultra-F Light-Activated, Longer-Term Temporary Material

- Perfect for longer term temporaries.
- Rigid, harder wearing material, with fluoride release.
- High compressive strength for bruxism patients.
- Ideal for tooth-colored restorations in deciduous teeth.
- Less soluble than reinforced ZOE intermediate restorative material.
- Light-activated, sets in 40 seconds or less.
- No etching or bonding agent required.
- **Tempit Ultra: without fluoride**



NoMIX® Short-term, Temporary Cement

- Easy to use, dispense, clean and remove from prep.
- Low film thickness, harder set.
- Tube applicators are ideal for both chair-side application and for patient take-home.
- Expands slightly upon setting for a tight margin seal.
- Moisture activated, no need to dry.
- **Also available in syringe delivery for chair-side application.**



Access Crown® Temporary Crown & Bridge Material


- **More:** shades (7), speed, strength and savings.
- **Ideal handling characteristics:** ready to remove from mouth in 60 seconds.
- Low shrinkage, low heat generation.
- High compressive, flexural strength.
- Highly aesthetic – polishes easily, beautifully.




Natural Esthetics, High Fluorescence, Polishes Easily.

✓ MORE SPEED

Intra Oral Set & Final Cure Time



✓ MORE STRENGTH



PARTNER WITH US!
CENTRIXDENTAL.COM
800.235.5862



MORE OPPORTUNITIES FOR
PATIENT & PRACTICE HEALTH!



LEARN HOW WE MAKE
DENTISTRY EASIER
FOR YOU!

DELIVERY SYSTEMS – CONTROL COSTS AND ELIMINATE CROSS-CONTAMINATION



Benda® Brush & Benda® Micro Disposable Brush & Micro Applicators

- The original Disposable Bendable Applicator.
- Benda Micro – use when you need to apply only one drop.
- Benda Brush – use when you need to apply more than a drop.
- Regular point for large applications.
- Mini point for smaller spaces.
- Made in the USA.



Snap Fit™ & Snub Nose Snap-Fit™ Composite Resin Guns

- Easy, direct-to-prep delivery of today's dental materials.
- Easy loading.
- Exclusive handle design.
- Works with all tubes and plugs.
- Made in the USA and used in every dental office!



AccuDose® & NeedleTube® Delivery Systems

- Dispense exactly the material needed.
- Prevent cross-contamination.
- Ready to use, dispense directly into the prep.
- No waste, no contamination.
- No more 'composite keeper'!
- Compatible with all C-R Syringes and guns.
- Made in the USA.

POST-ENDO TREATMENT – DUAL-CURE BUILDUPS THAT BOND



Connexio® Dual-Cure, Self-Etch Bond for Core Buildups


- Designed for the demands of deep core build-ups where light cannot reach.
- Dual-cure adhesive system provides the highest possible bond strength.*
- Nanoparticle reinforced.
- Water solvent (not acetone or ethyl alcohol); makes Connexio more tolerant of humid canal conditions than other adhesives.

Encore® D/C MiniMix® Dual-Cure, Core Buildup Composite

- Non-slumping formula stays where it's placed, will not run and does not require a matrix or core form.
- Exclusive MiniMix® cartridge (0.5 mL) controls costs while eliminating waste and cross-contamination.
- Cuts like dentin to reduce ditching.
- NeedleTube® mixing tip delivers material to bottom of the post hole.
- Made in the USA.



The First Single-Patient Cartridge!



IMPRESSIONING – EASIER, FASTER BITE REGISTRATION

Access BLUE® VPS Quick Bite Registration Material

- Mousse-like viscosity for easy, fast application.
- Sets hard, like laboratory plaster in just 30 seconds.
- Most accurate bite registration for intercuspal and centric relationships.
- Cut and trim the impression without it crumbling.
- Made in the USA.

EARN CE CREDITS WITH FREE LUNCH & LEARNS!

- No More "Watch & Wait" Treat & Restore with Silver Diamine Fluoride
- Prevention for Life – Getting Your Patients Through Life Caries Free
- The Core Issue – Structural Integrity of Broken Down Teeth
- How to Avoid Common Mistakes When Taking Impressions
- A Simplified Approach to Temporary Crown & Bridge

To schedule a Centrix Lunch & Learn, call 800.235.5862 or visit centrixvirtuallearning.com/lunch-and-learn.

CENTRIX INCORPORATED

770 River Road | Shelton, CT | 06484 | USA
203.929.5582, Toll-Free (USA & Canada): 800.235.5862

Centrix is an ISO 13485 certified dental manufacturer and our products are marketed worldwide with FDA, CE, Health Canada, TGA and PAL registrations on file, 0284. © 2024 Centrix, Inc. | MK0203 RevB

