

A DIFFERENT VPS. MAKING IMPRESSIONING EASIER!

HERE'S WHAT MAKES NOCORD DIFFERENT AND BETTER:

- NoCord VPS is the first VPS with built-in hemostasis.**
Stops bleeding and fluid flow for clear margins, exceptional detail.
- NoCord VPS MegaBody™ Tray Material creates sulcular expansion.**
Drives the NoCord Wash Material into the sulcus, displacing soft tissue, exposing all margins.
- NoCord VPS is rigid and stable.**
Great results in the most demanding cases.
- NoCord VPS compliments your technique.**
Whether you use two cords, paste or no retraction at all, NoCord makes impressioning safer, easier and more error free.



Clean Capture. Beautiful Results!


HERE'S WHAT MAKES NOCORD BETTER FOR YOU:

OFTEN: YOU WILL NEED NO CORD.




Supragingival margins and relatively healthy tissue is the perfect case for NoCord!

MANY TIMES: ONLY ONE CORD.



The sulcular expansion NoCord creates often eliminates the need to place a second cord.

EVERY TIME: SAFER, EASIER, BETTER.




NoCord Wash and its 15% alum (ammonium aluminum sulfate) astringent against the gingiva helps eliminate re-bleeders and pooling of blood, while protecting fragile tissue.

NOCORD™ VPS MAKES IMPRESSIONING:

SAFER.

- Gently expands the sulcus.
- No cord or packing instrument.
- Protects against damage to the epithelial attachment.
- No concerns about recession.



“When a patient has **thin** tissue that is **prone** to recession or recently had a periodontal procedure, I can **avoid packing cord** and safely record the detail.”


David Polhemus, DDS, AACD

“...easily take highly accurate impressions while **causing no harm to sensitive gingival tissues**. NoCord VPS affords a clean capture without extra tissue retraction.”

Summer Lerch, DDS, FAGD, AAACD

EASIER.

- NoCord Wash. Integrated hemostasis.
- MegaBody 50 mL cartridges that are easier to extrude.
- 380 mL cartridges for automatic mixing machines.
- Dual arch or full arch, your preference.



“NoCord is a trifecta for impression materials; **sulcus expansion, hemostasis and accuracy**. NoCord is the Next Generation of impression materials.”


Frank Milnar, DDS



WORK ≤1:00 min SET ≥4:45 min

BETTER!

- Single units or large span cases, NoCord provides an extra margin of safety.
- Rigid. More accurate implant impressions.
- For cord users, controls rebleeders when a cord is removed.
- Scans easily for detailed digitization.



“**The stiffness of the materials** provides an accurate transfer of the implant impression coping while also **recording excellent soft tissue detail**.”

David Polhemus, DDS, AACD

“The use of NoCord impression material, combined with its **hemostatic agents**, can also **stop the oozing** and safely widen the **sulcus** to allow a more accurate and detailed capture of the margins. Having this option can help **prevent the need to have the patient return** for another appointment to attempt a better intraoral scan or more time consuming use of cords or pastes at that appointment.”

Gregori Kurtzman, DDS, & Tarun Agarwal, DDS

ORDERING INFORMATION

- REF 310280** NoCord VPS Trial Kit
(1 x 50 mL Wash Cartridge, 2 x 50 mL MegaBody™ Tray Material Cartridge, 24 x SuperMixer™ Mixing Nozzle)
- REF 310281** NoCord VPS MegaBody Tray Material 50 mL, 2-pk.
(2 x 50 mL Cartridge, 6 x SuperMixer Mixing Nozzle)
- REF 310282** NoCord VPS Wash Material 50 mL, 2-pk.
(2 x Cartridge 50 mL, 6 x SuperMixer Mixing Nozzle)
- REF 310283** NoCord VPS MegaBody Tray Material 50 mL, 10-pk.
(10 x 50 mL Cartridge, 12 x SuperMixer Mixing Nozzle)
- REF 310284** NoCord VPS Wash Material 50 mL, 10-pk.
(10 x 50 mL Cartridge, 12 x SuperMixer Mixing Nozzle)

- REF 310285** NoCord VPS System Intro Kit
(6 x 50 mL Wash Cartridge, 12 x 50 mL MegaBody Tray Material Cartridge, 24 x SuperMixer Mixing Nozzle)
- REF 310286** NoCord VPS MegaBody Tray Material 380 mL
(1 x 380 mL Cartridge, 10 x SuperMixer, 1 x Bayonet Ring)
- REF 310287** NoCord VPS MegaBody Tray Material 380 mL, 2-pk.
(2 x 380 mL Cartridge, 1 x Bayonet Ring)
- REF 360005** Access™ Tray Adhesive, 10 mL bottle
- REF 360010** SuperMixer Mixing Nozzles, Teal Hub, 24
- REF 910631** SuperMixer Mixing Nozzles, Yellow Hub with 18 ga Green Tip, 50
- REF 910637** SuperMixer Mixing Nozzles, 380, Yellow Hub, 40
- REF 910638** Bayonet Ring, Yellow, 2

Ordering is easy! Call **800.235.5862** or contact your preferred dealer. Visit **centrixdental.com** for more information

Photo credits: #1, 2, 4 & 8; Dr. Miguel Ortiz, Brookline, MA. #3; Dr. Summer Lerch, New Haven, CT #5 & 6; Dr. Amit Desai, Cheshire, CT. #7, 9, 10; Alternative Dental Lab (Dental Wings Software), New Haven, CT. #11; Dr. Cristian Enachesau, Romania.

© 2023 Centrix, Inc. All rights reserved. NoCord VPS, NoCord Wash Material, NoCord MegaBody Tray Material, and SuperMixer are trademarks of Centrix, Inc. | MK0099 REV F

centrix
MAKING DENTISTRY EASIER®

MAKING DENTISTRY EASIER® | CENTRIXNOCORDVPS.COM